





Glunovo E2 APP

User manual

Compatible with Glunovo i3 & p3 CGMS

Agenda

- #1 What's new?
- #2 Compatible Smart Phones
- #3 Downloads Applications
- #4 App Usage Tips
- #5 App Use Guide
- #6 Glucose Data Uploading



	Glunovo E1	Glunovo E2
Compatible Hardware	i3	i3 & p3
Launch Time	2021.2	2022.11
Latest Version	V1.3.13	V1.1.1
Operating Strategy	Stop Maintenance at 2023 Q4	2022.11
Compatibility with other software	Support GN-Share, GN-Carer	GN-Share Ver.1.4.0 GN-Carer :2022 Q4



If you are using glunovo E2, you must use GN-Share Ver.1.4.0 for compatibility

#1 What's new ? –App Icon

E1



Glnovo
Utilities

OPEN



GN-Share App V1.1.0

- Get data from Glnovo App
- Sync to the cloud
- Manage the observer list

E2



GlnovoE2
Utilities

GET



GN-Share App V1.4.0

- Get data from Glnovo App
- Sync to the cloud
- Manage the observer list

#1 What's new ? –p3 package



Brand	Model	App Ver.
Huawei	mate 20	10
Redmi	K30Pro	10
Samsung	SM-S8	9
	SM-A53	12
	Galaxy S21 FE	12
Real me	GT 2	12
Apple	iPhone X	15.5
	iPhone 11	15.6
	iPhone SE	14.4


If your smartphone is not listed, it does not mean you cannot use your device with Glunovo product. We are continuously testing more devices for compatibility



#2



Smart Phones under testing

Brand	Model	System
Apple	iPhone 7	IOS 14.2
Apple	iPhone XS Max	IOS 14.4.2
Apple	iPhone 8	IOS 14.6
Honor	HRY-AL00a	Android 10
Red MI	Redmi 9A	Android 10
Red MI	Redmi 10 pro	Android 11
Samsung	Galaxy A70	Android 9
XIAOMI	XIAO MI10 pro	Android 11
XIAOMI	XIAO MI9	Android 11






#3 Download the app



Transmitter	Sensor Started	Last
CAL0036	2022-11-18 15:31	2022-11-18 15:31
E5M0246	2022-09-10 00:19	2022-09-10 00:19
I5M0263	2022-09-10 00:00	2022-09-10 00:00



#4

App Configuration Tips

We recommend users configure smart phones as shown to keep your Glunovo App functioning properly

<https://www.infinovo.com/en/lists/58.html>




#4 App Configuration

iOS Device Setting Path

Setting -> Glunovo E2

iOS Configuration Tips

- 01-Allow Bluetooth Access
- 02-Allow Notifications
- 03-Enable Background App Refresh





01-Turning-off Low Battery Mode

Path: Setting -> APPS-> Glunovo E2-> Battery-> Unrestricted

Change Setting to **Unrestricted**

02-Enable App Background Refresh

Path: Setting -> General Management->Battery->Background usage limit->Auto run apps

Enable Glunovo E2 to auto run apps

03- Enable Push Notification


Path: Setting -> Notification->Glunovo E2

Allow Glunovo E2 notifications

04- Do not Add Glunovo to deep sleeping apps


Path: Setting -> General Management->Battery->Background usage limit->Deep sleeping apps

DO NOT ADD Glunovo E2 TO Deep sleeping apps



#5


Glunovo App User Guide



#5 App User Guide

01- Fill the personal information form & Agree 《Use Terms and Privacy Policy》

02- Tab on **[Next Step]**
Start using Glunovo E2



Glunovo[®]

Personal Information


First Name

Last Name

Country

Next Step

Agree 《User Treaty and Privacy Policy》



17:28

Glunovo

INFINOVO TERMS OF USE

THESE INFINOVO TERMS OF USE ARE EFFECTIVE AS OF November 20, 2020.

These Terms of Use constitute a legally binding agreement made between Infinovo ("Infinovo", "we", "us", "our") and you, whether personally or an entity for whose use and on whose behalf you purchase a Infinovo Product or service ("you") , concerning your use of (including any access to): the Infinovo website located at www.infinovo.com and other websites we operate (our "Infinovo Websites"); the software applications we make available for download or access at our Websites, in the Google Play or App Store or other mobile service provider locations or other locations we indicate ("Apps"); and any Infinovo products ("Infinovo Products") and our data services available through the internet to users of Infinovo Products or Software Apps that provide and permit the access, collection, storage, processing, analysis and/or transmission of data generated by a Infinovo Product or Software App ("Infinovo Services"); as well as any email or other messages that you receive from us. **YOU ARE ENTITLED TO USING INFINOVO WEBSITES, INFINOVO PRODUCTS AND INFINOVO SERVICES ONLY AFTER YOU AGREE WITH THESE TERMS OF USE.**

All of the following terms are a single legal agreement between you and Infinovo.

The information provided on the Infinovo Websites is not intended for distribution to or use by any person or entity in any jurisdiction or country where such distribution or use would be contrary to law or regulation or which would subject us to any registration requirement within such jurisdiction or country. Accordingly, those persons who choose to access Infinovo Websites from other locations do so on their own initiative and are solely responsible for compliance with local laws, if and to the extent local laws are applicable.

Infinovo Websites, Products and Services are not for medical decisions.

You understand that Infinovo is not a healthcare professional and does not provide medical, health or other professional services or advice, nor do we verify the authentication and accuracy of User Data. Infinovo Products and Services are not replacement for proper medical care, and you agree that the user is solely responsible for obtaining proper treatment for his or her conditions. Certain information on the Websites and Apps relates to products and services offered by Infinovo, which is for your general education and information only. Information you read on the Websites and Apps cannot replace the relationship between a health care professional (e.g., a physician) and a patient. If you think that you are having a medical emergency while using our product, refer to your doctors and healthcare professionals before you make any medical decisions.

- Do NOT use or rely on Infinovo Product and Services or any associated materials for any medical purpose or decision.
- Do NOT rely on Infinovo Product and Services for any real-time alarms or time critical data.
- Do NOT use or rely on Infinovo Product and Services for treatment decisions or use as a substitute for professional healthcare judgement.
- Any part of the system can fail at any time. Always seek the advice of a qualified healthcare professional for any medical questions.

You may provide the information and reports received from Infinovo Product and Services to certain third parties selected by you at your own responsibility.

Intellectual property rights

Unless otherwise indicated, the Infinovo Products, Services and Websites are our proprietary property and all resources code, databases, functionality, software, website designs, audio,

First-time user of Glunovo E2 has to set information below:

01-Set Unit of Glucose

The user needs to select the glucose unit. All glucose Settings and information in the APP are displayed according to the selected unit. Once the glucose unit is selected, it cannot be modified ;

02-Set hyperglycemia threshold & Alert

User needs to set the hyperglycemia and alert thresholds, and set whether to enable the hyperglycemia alert and the alert mode.

03-Set hypoglycemia threshold & Alert

User needs to set the hypoglycemia and alert thresholds, and set whether to enable the hypoglycemia alert and the alert mode.

04-Set threshold of urgent hypoglycemia

International uniform threshold for emergency hypoglycemia, which cannot be modified by users; When the user's glucose level is lower than the threshold, it means that the user's life will be threatened, and the user will be informed by the strongest reminder.



#5 App User Guide -Transmitter Pairing

At transmitter pairing screen, fill out Transmitter SN and Pin on the Label;

Tab **START PAIRING**, Start pairing the transmitter;

After successful pairing, navigate to sensor pairing screen;

Transmitter SN
O5M0014 - BED5

Transmitter pin
BED5

Glunovo[®]

PLEASE ENTER TRANSMITTER SN

The Transmitter SN is stamped on the Transmitter and also on the CGM packaging

O5M0014

TRANSMITTER PIN

The Transmitter PIN is stamped on the Transmitter and also on the CGM packaging

BED5

START PAIRING

Need Help?

Transmitters of Continuous Glucose Monitoring Systems

REF	CP04020201	#	p3	Qty	1pc
	XXX-X-XX-XX		XXXX-XX-XX		
SN	O5M0014	LOT	XXXXXXXXXX		

UDI

(01) 0 6973 18 10018 2
(10) XXXX XXXX
(17) 000000
(21) XXXXXXXX-XXXX

Infinovo
Innovation for Diabetes


Infinovo Medical Co.,Ltd.
3rd Floor,6th Building
No.888,Zhujiang Road
Rudong
226400 Jiangsu
China

EC REP Lins Service & Consulting GmbH
Obere Seegasse 34/2,69124
Heidelberg, Germany
TEL: 0049 175 4870 819

32°F (0°C) 113°F (45°C) 10% 95%

MD CE 0197

! We do not recommend to skip factory calibration code




LOT: CP04010201

Sensor Code: 0720



P3 Sensor (New)




I3 Sensor

After the device is connected to the transmitter, the system automatically determines whether the device belongs to i3 or p3 and displays different screen based on the device model. If the p3 device is connected, the user can choose to skip the factory calibration code, but calibration is required after startup ;


- Start-up time:
- I3: 190mins
 - P3: 120mins

! When glucose reading does not match the body's perception, please take fingerstick







The image shows a smartphone displaying the app's home screen. The screen features a large glucose reading of 2.7 mmol/L, a trend graph for the last 4 hours, and a bottom navigation bar with icons for Home, Calibration, Report, Logbook, and Setting. Callout boxes with dashed orange borders point to specific elements on the screen, providing detailed descriptions of each.

- Latest glucose data**: Points to the large '2.7' reading.
- Time of getting the latest data**: Points to '3 mins ago'.
- Device expiration time**: Points to '5 days left'.
- Glucose trend graph**: Points to the green line graph, with a sub-description: 'Display user's glucose trend in 4 hours'.
- Tab menu**: Points to the bottom navigation bar.
- Add a new event**: Points to the '+' icon in the top right.
- Glucose variation trend**: Points to the arrow icon next to the glucose unit.
- Glucose unit**: Points to 'mmol/L'.
- Display events in 4 hours**: Points to the 'Suspend 15 mins' notification bar.



Glucose variation trend Icons and definition


Icon	Definition	
→	Steady	Glucose does not rise or decrease by more than 0.06 mmol/L (1.08 mg/dl) each minute
↗	Rising Slowly	Glucose could rise 0.06-0.11 mmol/L (1-2 mg/dl) each minute
↑	Rising	Glucose could rise 0.11-0.17 mmol/L (2 -3 mg/dl) each minute
↑↑	Rising Rapidly	Glucose could rise more than 0.17 mmol/L (3 mg/dl) each minute
↘	Falling Slowly	Glucose could fall 0.06-0.11 mmol/L (1 -2 mg/dl) each minute
↓	Falling	Glucose could fall 0.11-0.17 mmol/L (2 -3 mg/dl) each minute
↓↓	Falling Rapidly	Glucose could fall more than 0.17 mmol/L (3 mg/dl) each minute



Users can record the events in daily life that will affect the change of glucose, such as insulin, meals, drugs, exercise, etc., as well as the finger blood;


After the event is saved, the corresponding event will be displayed in the diagram on home screen.

#5 App User Guide – Home screen (cgm information)




Tab **CGM** to see detailed device information:

- Operating status of equipment
- Serial number of transmitter
- Batch number of the sensor
- Sensor current
- Last calibration event
- App Version



p3 is calibration-free CGM. If the deviation between the product and the measurement result of fingerprick is too large, user could choose to calibrate. Please refer to user guide for details :

- First and second days of wearing, deviation tends to be larger than the rest of the days, user do not need to calibrate;
- If you need to calibrate, please calibrate when glucose is steady ;
- After taking fingerprick, please input the BG reading immediately for calibration;
- We do not recommend users calibrate more than 2 times/day;



The report displays user' s glucose analysis report of current round.
The analysis report has two forms:

- **TIR: Time In Range**
Demonstrates distribution of different glucose levels in target time range
- **Glucose Profile**
Displays glucose distribution in 24 hours




#5 App User Guide – Logbook



Filter:


- Categories: Sport, Medication, Carbs, Insulin;
- p3 CGM reading;
- Hyperglycemia/hypoglycemia Alert;
- Calibration Record;
- BG Record;

Swipe left to edit or delete an event or BG record:



- View device information (Sensor an transmitter)
- User can remove current device by selecting **【Remove Current Device】**

App version & Copyright





Export historical glucose data and event

Set/Edit hyperglycemia and hypoglycemia threshold and alert

#6

Glucose Data Uploading









Remark:

1. When uploading data from Glunovo E2 APP, please **update GN-Share APP to V1.4.0 version**, otherwise user will not be able to upload blood glucose data from Glunovo E2.
2. Please **keep Glunovo E2 and GN-Share APP running in the background** in order to ensure user timely uploading of glucose data.
3. GN-Share APP can be connected to Glunovo E1 and Glunovo E2 at the same time, but **only one APP can be used to collect data and upload at the same time**; for regions where i3 and p3 are sold at the same time, please instruct users to uninstall Glunovo E1 and use Glunovo E2.

#6

Glucose Data Uploading



Please keep the apps running in the background.

Continuous improving strategy

➤ Improving Glunovo E2

- Improving Bluetooth Connection
- Keep testing latest smart phones for compatibility
- Improving App's ability of keeping active

